PHYSICIAN INSTRUCTIONS

GRAETHER 2000™
PUPIL EXPANDER SYSTEM

STERILE R

Rx ONLY

 **EagleVision**®
A KATENA BRAND

Description

The Graether 2000 Pupil Expander System, REF1550, consists of three components: Pupil Expander, Iris Glide-Retractor, and Insertion Instrument. (figs. 1-3)

The **Pupil Expander** (fig. 1) is a clear silicone ring with an external diameter of 7.5mm and an aperture of 6.3mm. The ring has a gap of approximately 3.5mm in its circumference that is bridged by a strap. The strap terminates at both ends into tabs with holes on the upper surface of the ring. The ring has a groove or slot in its outer circumference (fig. 1) that accepts the iris sphincter.

The **Iris Glide-Retractor** (fig. 2) is plastic, self-retaining and positioned in the primary incision prior to placement of the Pupil Expander ring. It serves to retract and fix the sphincter at the internal margin of the incision, providing counter-traction as the Pupil Expander is inserted.

The Pupil Expander is supplied ready to insert on a pre-loaded **Insertion Instrument** (fig. 3) with the Pupil Expander ring stretched around the tip. The instrument allows insertion of the Pupil Expander through a 2.75mm* or larger incision, gliding on the surface of the Glide-Retractor. When positioned within the stretched pupil, the Pupil Expander is released by pressing the trigger button on the insertion instrument.

***Note:** If an incision smaller than 2.75mm is used, the tab on the Glide-Retractor can be cut off with scissors to allow the Glide-Retractor to pass into the anterior chamber.

Indications For Use

The Graether 2000 Pupil Expander System is designed to provide an enlarged pupil opening during phacoemulsification and lens implantation. It expands the pupil to 6.3mm and sustains and protects the pupil margin during the operation. Its soft silicone wall applies uniform pressure against the sphincter and is atraumatic to delicate ocular tissues. Over-stretching and point traction on the sphincter is avoided, minimizing the risk of sphincter paralysis. The gap in the ring gives access to the pupil aperture while the strap protects the iris from chafing. The ring maintains a constant relationship to the capsule opening, facilitating the capsulotomy and accurate placement of the lens haptics.

Contraindications

The Graether 2000 Pupil Expander System is contraindicated for patients with known silicone or plastic sensitivity.

Capsulotomy and Hydrodissection

A continuous circular capsulotomy can be performed with a cystotome and/or forceps using the margin of the Expander as a template. It is recommended that a fairly large capsulotomy be made to facilitate cataract removal and lens implantation and to prevent pupil block. *If hydrodissection of the nucleus is employed, it is important that care is taken to avoid entrapment of fluid within the capsule. Injection should not be too rapid and the nucleus should be tilted during injection to allow the escape of fluid around it. On rare occasion, the uniform pressure of the Pupil Expander ring may impede the escape of fluid from the capsule and produce a form of pupil block. If such a block occurs, the chamber will shallow and the eye will become firm. It can be easily relieved by backward pressure on the nucleus to allow the irrigating fluid to escape.* After release of the fluid, the chamber should be re-deepened with visco-elastic prior to phacoemulsification.

Phacoemulsification and Lens Implantation

The surgeon's preferred technique can be employed. Tumbling of the entire nucleus should be avoided as the Pupil Expander may limit the forward movement of a large, firm nucleus and place undue stress on the posterior capsule. Placement of the lens implant is facilitated because the edge of the capsule can be easily identified.

NURSE

Check the contents immediately prior to surgery:

Peel open the Tyvek tray lid over the sterile field, keeping the tray level with the floor. (fig. 4)

The blue disc covering the small Glide-Retractor is removed with forceps to expose the Glide-Retractor in the tray well. It is suggested that the blue disc be placed in the staging area and the Glide-retractor be removed with a smooth forceps and placed on top of it until needed.

The inserter is staged in the tray to prevent premature release of the Pupil Expander.

Preparation

All pre-operative and intra-operative pharmacologic aids, including topical antiprostaglandins, cycloplegics and midriatics are employed prior to surgery. Adrenaline may be added to the irrigating solution used during surgery, but is not essential. Retrobulbar or peribulbar anesthesia is not necessary. If topical anesthesia is used, intracameral lidocaine without preservative is recommended. Either a scleral or clear corneal incision is satisfactory. Any visco-elastic material may be employed. If the chamber is shallow or if there is corneal disease, a more viscous substance is recommended. Pre-stretching of the pupil is not necessary, but lysis of dense synechiae or removal of any inflammatory membrane from the pupil should be done prior to the placement of the Pupil Expander.

The Glide-Retractor is removed from the staging area with a smooth forceps. It may be helpful to grasp the Glide-Retractor with forceps jaws aligned with its axis to enhance stiffness for ease of insertion. The incision is opened with a toothed forceps or thin, flat spatula like the Jaffe spatula. The Glide-Retractor tip is passed through the incision and inserted so that the tab on the right side passes into the anterior chamber. In a very small incision, the spatula will help guide the tab past the incision by sliding the Glide-Retractor *under* the spatula. (fig. 5) The hook at the tip of the Glide-Retractor engages the rim of the iris sphincter and the Glide-Retractor is withdrawn until the tab catches the edge of the external incision. (fig. 6) The engaged sphincter will be pulled toward the incision. The anterior chamber is filled with visco-elastic.

The pre-loaded Pupil Expander on its inserter is removed from the sterile delivery tray.

Note: *If the Pupil Expander has been displaced from the instrument during shipment or handling, it can be reattached under the microscope. The instrument and expander are turned over exposing the undersurface. (The post of the instrument is positioned up.) Place the Pupil Expander ring on a flat, sterile surface with the tabs against the surface. Using a smooth forceps, approach the strap from inside the ring. Grasp the strap at its midpoint. Center the opposite end of the ring on the tip of the instrument and loop the strap over the metal post, making sure the undersurface of the strap is toward the post and the tabs are aligned. (fig. 7)*

INSERTION

Overview

There are six easy steps to insert the Pupil Expander:

1. Insert the Glide-Retractor and retract the pupil toward the incision.
2. Insert the Pupil Expander and push the pupil margin to the opposite limbus while retracting the Glide-Retractor.
3. Make sure the pupil is fully stretched by aligning the tabs of the Pupil Expander just past the folded end of the Glide-Retractor.
4. Push the button to release the Pupil Expander.
5. While holding the Pupil Expander ring in place with the tip of the insertion instrument, withdraw the Glide-Retractor from the eye.
6. Center the pupil and withdraw the insertion instrument from the eye.

The incision is opened with toothed forceps, and the insertion instrument is inserted through the incision, traveling on the surface of the Glide-Retractor. (fig. 8) The inserter tip engages the sphincter opposite the Glide-Retractor and the pupil margin is pushed only as far as the opposite limbus. (fig. 9) Do not force the inserter tip beyond the limbus. The toothed forceps is used to simultaneously retract the sphincter toward the incision by traction on the Glide-Retractor. (fig. 9) This counter traction stretches the sphincter between the inserter tip and the Glide-Retractor enough to accommodate the Pupil Expander ring.

When the tabs on the Expander ring are positioned beyond the Glide-Retractor tip (fig. 10), the Pupil Expander is released by pressing the trigger button on the handle. The Pupil Expander ring will start to expand onto the iris sphincter. (fig. 11) While holding the Pupil Expander ring in place with the tip of the insertion instrument, withdraw the Glide-Retractor from the incision. (fig. 12) (Some resistance will be felt as the tip unfolds against the incision.) The released iris sphincter will move into the slot in the Pupil Expander ring as the ring resumes its natural circular shape. The Expander ring will cover the edge of the sphincter, with the tabs resting against the iris and the strap will bridge the sphincter adjacent to the incision. (fig. 13) The pupil can be centered with the tip of the insertion instrument prior to withdrawing the insertion instrument from the eye. (See warnings.) If a portion of the Expander ring fails to engage the sphincter, it can be easily nudged into place with the tip of the spatula positioned in the slot of the ring or with the tip of the visco-elastic canula used to refill the chamber. (fig. 14)

REMOVAL

Removal of the Expander ring is simple. The strap is pushed into the pupil to disengage the ends of the ring from the iris. (fig. 15) The strap is grasped and the ring withdrawn. (fig. 16) Alternatively, the strap can be cut with scissors and one end of the ring lifted from the iris and withdrawn.

Conditions Encountered and Their Solutions

- **Difficulty engaging the iris sphincter.**
The sphincter must be fully stretched before releasing the Expander from the insertion instrument. Make sure the tabs on the Pupil Expander ring have cleared the folded tip of the Glide-Retractor, before releasing the Expander. Hold the insertion instrument in place as the Glide-Retractor is withdrawn from the incision. Resistance will be encountered as the tip of the Glide-Retractor unfolds against the incision.
- **Occasionally one end of the ring will not be fully engaged or the strap will remain partially folded.**
The free end of the ring can be pushed into the pupil space with the tip of a Jaffe type spatula and allowed to expand onto the sphincter. Rarely two instruments are necessary to further stretch the ring or strap. (fig. 17)
- **Shallowing of the anterior chamber or a sudden increase in IOP.**
This is almost always due to a relative pupil block as fluid is injected for hydrodissection. Release fluid by pressing back or tilting the nucleus.
- **Phaco tip or sleeve displaces the Expander ring from the sphincter.**
This can be avoided by making sure the strap is positioned toward the iris even if not fully stretched. A small amount of visco-elastic can be placed above the strap to expand the space just before inserting the phaco tip. Alternatively, the strap can be cut or removed.
- **A bound down sphincter or rigid pupil.**
If the pupil will not dilate at all with miotics, it may be possible to strip a membranous scar from the sphincter muscle. Free up the pupil with a collar button retractor or sideport manipulator. Grasp the membrane with a Utrata or similar capsulotomy forceps and gently pull it tangentially from the sphincter. It is not necessary to pre-stretch the pupil prior to inserting the Pupil Expander. Occasionally a pupil will be encountered that is so rigid that the sphincter must be cut to avoid tearing it. A sphincterotomy can be performed and the Pupil Expander applied to dilate the pupil and control the leaves of the iris. If necessary, the sphincter defect can be closed with a suture.
- **A nucleus too firm to emulsify.**
These exceptional cases are best handled by the surgeon's standard extracapsular technique. The Pupil Expander can be used to facilitate a large continuous circular capsulotomy, but it should be removed before expressing the lens nucleus. Often a sphincterotomy is desirable in this situation rather than attempting to preserve the sphincter.

Cautions

U.S. Federal law restricts this device for sale by or on the order of a licensed medical practitioner.

Warnings

The Graether Pupil Expander System is specifically designed as a safe, efficient method for expanding a pupil that will not otherwise dilate. Using both the Insertion Instrument and the Iris Glide Retractor, as described in the Insertion Instructions, ensures a successful result.

If the operating Physician determines that an alternative insertion technique is to be used, one omitting the use of the IGR, the following precaution should be observed. (fig. 18)

The inserter instrument incorporates a pin protruding from the ventral rear of the tip which holds the expander ring in place. Be aware that the pin may contact and cause injury to the tissue. Pressing the release button withdraws the pin into the inserter tip releasing the expander ring. After release of the expander ring, some portion of the pin may be left exposed. DO NOT allow the pin to PRESS against the tissue before, or during, insertion of the expander ring.


















- Single Use Only. Katena Products, Inc. assumes no liability for products which have been resterilized by health care facilities. The device may be compromised if it is reused due to damage when explanted or degradation due to resterilization. There is also a potential of infection.
- DO NOT RESTERILIZE.
- Supplied Sterile. Sterility not guaranteed if package has been opened or damaged.

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