



blink SINGLE-USE
INSTRUMENTS

Single-use Ophthalmic Instruments

OPHTHALMOLOGY

 **corzamedical**



For ophthalmic surgeons and healthcare facilities seeking the highest standards in patient safety and surgical precision Corza Medical provides an extensive range of ophthalmic single-use instruments.

Unlike traditional reusable instruments, our Blink single-use instruments eliminate the risk of cross-contamination, ensuring sterility with every use. Blink single use ophthalmic instruments are meticulously designed for superior performance and convenience, providing a reliable, cost-effective alternative that enhances both patient outcomes and operational efficiency. Choose Corza Medical for innovative, trustworthy, and ready-to-use ophthalmic instruments that elevate the standard of care.

| | |
|------------------------------------|----|
| Tying Forceps | 4 |
| General Forceps | 4 |
| Fixation Forceps | 5 |
| Jewelers Forceps | 6 |
| Mosquito Forceps | 6 |
| Capsulorhexis Forceps | 7 |
| Sprung Scissors | 8 |
| Iris Scissors | 8 |
| Tenotomy Scissors | 8 |
| Speculums | 9 |
| Manipulators | 10 |
| Choppers | 10 |
| Double-ended Manipulators/Choppers | 11 |
| Needle Holders | 11 |
| Refractive | 12 |
| Corneal | 12 |
| Markers | 12 |
| Sub Tenon Anesthesia | 13 |
| Oculoplastics | 14 |

Manufacturing and quality

All Blink single-use instruments are manufactured to the highest standards by ISO 13485:2016 accredited manufacturing facilities. Following the initial stages of manufacturing, all products are subject to strict inspection criteria at our FDA registered facility in Solihull, UK to ensure they meet the exacting standards expected by our customers. Following quality control, all products are packaged in our Class 8 cleanroom prior to sterilisation and global distribution.



Customer service and distribution

At our Parsippany, NJ office, we are dedicated to delivering exceptional customer service. Our knowledgeable team provides personalized support and solutions, ensuring every interaction is met with professionalism and care. We prioritize understanding and addressing the unique needs of our healthcare partners, fostering trust and long-term relationships. At Corza Medical, your satisfaction and success in improving patient outcomes are our highest priorities, and we are unwavering in our mission to enhance healthcare through superior service and cutting-edge medical device solutions.



Forceps

Tying Forceps

Tying forceps are available in five different designs, providing clinicians with a choice of straight and angled tips in both standard and short handle designs. The tying platforms are precisely manufactured to provide maximum grip and have smooth edges to prevent damage to sutures.

- **Tying Forceps Straight**
- **10-9009** | 5mm Tying Platform, O/L 90mm
- **10-9008** | 6mm Tying Platform, O/L 110mm

- **Tying Forceps Angled**
- **11-6016** | 5mm Platform Angled, O/L 90mm
- **11-6015** | 5mm Platform Angled, O/L 105mm

- **Tying Forceps Curved**
- **11-6012** | 5mm Platform Angled, O/L 108mm



General Forceps



- **Adson Forceps Fine Toothed Tip**
- **10-1305** | Fine Toothed Tip, O/L 120mm

- **Adson Brown Forceps**
- **10-1307** | O/L 120mm

- **Max Fine Forceps Curved**
- **11-6009** | Smooth Pointed Tip, O/L 110mm

- **Moorfields Forceps**
- **10-1301** | 16mm Serrated Platform, O/L 107mm

- **Moorfields Suture Utility Forceps**
- **10-1311** | 14mm Narrow Serrated Platform, O/L 107mm

- **Epilation Forceps**
- **10-1300** | Angled Tip, O/L 91mm

Fixation Forceps

Multiple handle designs are available with both notched and toothed tips. All tip designs are manufactured to provide a high level of grip for stabilising the eye during surgery whilst minimising trauma.

Hoskins Forceps Straight

10-8008 | 0.3mm Notched, O/L 90mm

10-8009 | 0.3mm Notched, O/L 108mm



Pierce Colibri Forceps

11-7009 | 0.3mm Notched, O/L 84mm



Bonn Forceps Straight

10-1106 | 0.25mm Toothed, O/L 106mm


10-1107 | 0.25mm Toothed, O/L 90mm



Colibri Forceps

11-7010 | 0.25mm Toothed, O/L 84mm

11-7011 | 0.25mm Toothed, O/L 108mm
(108mm version shown)



Lister Utility Forceps

10-1302 | 0.25mm Toothed, O/L 84mm

St Martins Forceps

10-1203 | 0.6mm Toothed, O/L 90mm

10-1204 | 0.3mm Toothed, O/L 90mm

Forceps

Jewelers Forceps

Available in five different designs, jewelers forceps provide clinicians with a choice for holding and manipulating fine tissues that can be easily damaged by other instruments.



Jewelers Forceps #1
10-2001



Jewelers Forceps #3
10-2003



Jewelers Forceps #5
10-2005



Jewelers Forceps #5 Angled
11-2005



Jewelers Forceps #7
11-2007

Mosquito Forceps

Mosquito forceps are available in both straight and curved options for grasping and mobilising small vessels during oculoplastic surgery.



Mosquito Forceps
11-1005 | Curved, O/L 100mm



Mosquito Forceps
10-1006 | Straight, O/L 100mm

Capsulorhexis Forceps




Capsulorhexis Forceps Utrata

- 11-5013** | Round Handle, Curved Shaft, Cystotome Tip, O/L 105mm
- 11-5012** | Round Handle, Straight Shaft, Cystotome Tip, O/L 105mm
- 11-5014** | Round Handle, Curved Shaft, Grasping Tip, O/L 105mm
- 11-5015** | Round Handle, Straight Shaft, Grasping Tip, O/L 105mm



Capsulorhexis Forceps Utrata

- 11-5050** | Flat Handle, Straight Shaft, Cystotome Tip, O/L 105mm
- 11-5051** | Flat Handle, Curved Shaft, Cystotome Tip, O/L 105mm
- 11-5052** | Flat Handle, Straight Shaft, Grasping Tip, O/L 105mm
- 11-5053** | Flat Handle, Curved Shaft, Grasping Tip, O/L 105mm



Cystotome Tip

A sharp tip design that can be used for creating the initial flap and grasping the lens capsule when carrying out a capsulorhexis.



Grasping Tip

Designed with larger grasping platforms for carrying out a capsulorhexis when the initial flap has been created with another instrument.


Capsulorhexis Forceps Inamura

- 11-5011** | Flat Handle, Straight Shaft, Cystotome Tip, O/L 105mm
- 11-5021** | Flat Handle, Curved Shaft, Cystotome Tip, O/L 105mm



Capsulorhexis Forceps Inamura

- 11-5020** | Round Handle, Curved Shaft, Cystotome Tip, O/L 105mm



Scissors | Speculums

Spring Scissors

Spring scissors are available in Westcott, Vannas and Ong designs. The Westcott designs have a variety of tip forms, as well as being available in straight or curved blade options. The drilled handle design increases grip and reduces the weight of the scissors, resulting in improved tactility.

Westcott Scissors

- 20-4013 | Straight, O/L 120mm
- 20-4015 | Straight, Sharp Tip, 120mm
- 21-4010 | Curved, O/L 120mm
- 21-4014 | Curved, Sharp Tip, O/L 120mm
- 21-4011 | Curved, Dolphin Nose, O/L 120mm
- 21-5001 | Curved, O/L 100mm



- **Vannas Scissors**
- 20-6000 | Straight, O/L 90mm
- 21-6001 | Curved, O/L 90mm



- **Ong Scissors**
- 21-3010 | Angled Tip
9mm, O/L 92mm


Iris Scissors and Tenotomy Scissors

- 20-7001 | Iris, Straight, Blunt Tip, O/L 100mm
- 20-7002 | Iris, Straight, Sharp Tip, O/L 100mm
- 21-7004 | Iris, Curved, Sharp Tip, O/L 100mm
- 21-7005 | Tenotomy, Curved, Pointed Tip, O/L 115mm
- 21-7006 | Tenotomy, Curved, Blunt Tip, O/L 115mm



Speculums

Speculums provide access to the eye for all types of ophthalmic procedures. The Blink speculum range includes a wide variety of designs. Most popular are Barraquer wire speculums, which are available in eight different designs and offer multiple benefits during ophthalmic surgery. Open Loop are widely used for cataract surgery, whereas Closed Loop and Solid Blade are favoured for administering ocular injections.



Barraquer Closed Loop

- 30-2004** | Straight, Width 30mm, Length 43mm
- 31-2003** | Temporal, Width 30mm, Length 42mm
- 31-2002** | Nasal, Width 30mm, Length 42mm



Kratz-Barraquer Open Loop

- 31-1006** | Temporal, Width 30mm, Length 40mm
- 30-1007** | Straight, Width 30mm, Length 42mm
- 31-1005** | Nasal, Width 30mm, Length 40mm
- 31-1009** | Temporal, Curved, Width 30mm, Length 42mm



Leibinger Speculum

31-3002 | Wire Loop, Temporal, O/L 65mm



Leibinger Speculum

31-3000 | Solid Blade, Temporal, O/L 65mm



Barraquer Speculum Solid Blade

30-4000 | Temporal, Width 30mm, Length O/L 44mm



Clarkes Speculum

31-3004 | Solid Blade O/L 95mm



Burratto Speculum

30-3001 | Wire Loop, O/L 40mm



Cook Speculum

31-3003 | Solid Blade, O/L 40mm



Paediatric Speculum Spring

30-3002 | Solid Blade Width 26mm, O/L 35mm



Neonatal Speculum Angled Temporal

31-3100 | Width 20mm O/L 35mm

Manipulators | Needle Holders



Manipulators

Available in ten different designs, the Blink manipulator range provides clinicians with a wide choice of options.



Sinsky Lens Hook Straight
40-5002



Sinsky Lens Hook Angled
41-5001



Nucleus Y Shaped Angled Horizontal Axis
41-2005



Nucleus Y Shaped Angled Vertical Axis
41-2006



Kuglen Hook Angled
41-4014



Mushroom Tipped Manipulator
41-8009



IOL Lester Manipulator
41-1007



Connor Wand
41-4013

Choppers

Available in eight different designs to meet clinicians' requirements when carrying out lens removal, these designs provide options for vertical and horizontal chopping techniques.



Drysdale Paddle Chopper
41-6007



Sadiq Chopper
41-4005



Koch Chopper 1.4mm
41-4008



Rosen Chopper 1.2mm
41-4010



Seibel Chopper
41-4012



Nagahara Chopper
41-4009

Double-ended Manipulators/Choppers

By bringing together our most popular designs in this range, these four options of double-ended manipulators/choppers provide the benefits of two products in one.



Sadiq Chopper/Mushroom Tipped Manipulator
41-9012



Sinsky Manipulator/Mushroom Tipped Manipulator
41-9010



Rosen Chopper/Mushroom Tipped Manipulator
41-9013



Nagahara Chopper/Connor Wand
41-9014

Needle Holders



- **Needle Holder Castroviejo**
- 90-7001 | 12mm Tying Platform, Straight With Lock, O/L 135mm
- **Needle Holder Castroviejo**
- 90-7002 | 12mm Tying Platform, Straight, O/L 135mm
- **Needle Holder Castroviejo**
- 91-7003 | 12mm Tying Platform, Delicate Curved, O/L 135mm
- **Needle Holder Troutman**
- 91-7004 | 6mm Tying Platform, Curved, O/L 105mm

Refractive | Corneal

Refractive



- **Toric Alignment Marker, 3, 6, 9 O'Clock**
- 91-6001



- **Toric Degree Gauge 11.7mm I.D**
- 91-6000



- **Toric IOL Marker 11mm Wide**
- 91-6002

Corneal

- **Epithelial Spatula**
- 51-0001

- **Epithelial Spatula, Rounded Tip**
- 51-0003

- **Alcohol Well**
- 51-0002 | 9.5mm

- **Stutman Spud and Bur**
- 90-1000



Markers

- **Double-ended Marker**
- 90-1015 | Flat Handle, 3.5mm-4mm, O/L 113mm

- **Braunstein Marker**
- 90-1017 | Round Handle, 3.5mm-4mm, O/L 80mm

- **Caliper & Rule**
- 90-1014 | 0-20mm Measurement Span



Sub Tenon Anesthesia

Corza Medical is focussed on providing healthcare professionals with innovative new products that save time, money and effort.

One such product is our Conjunctival Probe which is uniquely designed to create the perfect-sized incision for safe and effective injection of anesthetic during the Sub Tenon's procedure. The tip has been designed to be the same size as the Sub Tenon Cannula, resulting in far less backflow of anesthetic out of the wound, reduced red eye due to accidental puncturing of blood vessels, and improved akinesia due to the more effective distribution of anesthetic. The SISTA technique (Small Incision Sub Tenon Anesthesia) can be viewed below.


Should one of our standard sets not meet your requirements, please contact our helpful customer service staff who would be happy to help you specify your own custom pack.



• **Conjunctival Probe**
• CPO1



• **002 Sub Tenon Instrument Set**
• Scissors Westcott Curved Standard Tip
• Speculum Kratz-Barraquer Angled Temporal
• Moorfields Forceps



1 Raise the conjunctiva using a pair of non-toothed Moorfields forceps, holding them at 90 degrees to the direction of the probe.




2 Holding the probe parallel to the surface of the sclera, puncture the conjunctiva.




3 Keeping the probe in contact with the sclera, follow the curve of the globe for approximately 1cm.



4 Withdraw the probe but maintain the grip on the conjunctiva with the forceps.



5 Insert the Sub Tenon Cannula to inject the anesthetic in the usual manner.



• **270 CP Sub Tenon Set**
• Conjunctival Probe
• Speculum Kratz-Barraquer Angled Temporal
• Moorfields Forceps

Oculoplastics

Curettes

- **Curette Double-ended**
- **90-3007** | 2mm & 3mm Cups, O/L 130mm
- **Curette Double-ended**
- **90-3008** | 2mm & 5mm Cups, O/L 130mm
- **Curette Double-ended**
- **90-3009** | 2mm & 4mm Cups, O/L 130mm



- **Clamp Chalazion**
- **90-4005** | Small, 11mm Diameter Clamp, O/L 90mm
- **90-4006** | Medium, 18mm x 11mm Clamp, O/L 80mm
- **90-4007** | Large, 23mm x 15mm Clamp, O/L 90mm



- **092 – Chalazion Set – Small Clamp**
- **90-4005** | Chalazion Clamp Small
- **90-3009** | Curette Double Ended Small/Medium
- No.11 Scalpel
- Wrap 25x25cm

- **111 – Chalazion Set – Medium Clamp**
- **90-4006** | Chalazion Clamp Medium
- **90-3009** | Curette Double Ended Small/Medium
- No.11 Scalpel
- Wrap 25x25cm




- **Desmares Entropian Clamp**
- **90-4008** | Left, Clamp Diameter
- 24mm x 16mm, O/L 95mm



- **Desmares Entropian Clamp**
- **90-4009** | Right, Clamp Diameter
- 24mm x 16mm, O/L 95mm



- **Liebreich Probes**
- **90-2007** | Liebreich Probe 0000+000 | O/L 125mm
- **90-2008** | Liebreich Probe 000+00 | O/L 125mm
- **90-2009** | Liebreich Probe 00&0 | O/L 125mm
- **90-2010** | Liebreich Probe 1+2 | O/L 125mm
- **90-2011** | Liebreich Probe 3+4 | O/L 125mm



- **Thornton Ring**
- **90-1020** | 13mm Diameter



- **Desmarres Retractor**
- **90-9013** | 14mm Wide Blade, O/L 130mm



- **Graefe Strabismus Hook**
- **90-5004** | 10mm Tip Length, O/L 140mm



- **Alexandra Muscle Hook**
- **90-5005** | 11mm Tip Length, O/L 130mm



- **Backhaus Towel Clip**
- **90-5050**



- **Dilator Punctal Double-ended**
- **90-1011** | Pointed Tips, O/L 138mm




- **Dilator Nettle Ship**
- **90-1012** | 20mm Tapered Tip, O/L 118mm



- **Bulldog Clamp**
- **90-8012** | O/L 52mm

Standard Instrument Sets

Should one of our standard sets not meet your requirements, please contact our helpful customer service staff who would be happy to help you specify your own custom pack.



- **023: IVT Set**
- **90-1015** | Double Ended Marker
- **10-1104** | Forceps Suture Bonn 1:2 Teeth
- **30-1007** | Speculum Kratz-Barraquer



- **76-0073: Foreign Body Removal Set**
- **31-2003** | Speculum Barraquer Angled Temporal
- **90-1000** | Stutman Spud and Burr
- **10-2003** | Jewellers Forcep #3



- **210: Ophthalmic Suture Removal Set**
- **10-9007** | Forceps Mcpherson Tying Straight
- **21-6004** | Scissors Vannas Curved



- **76-0026A: Universal Cataract Set**
- **31-3002** | Leibinger Speculum
- **11-5050** | Forceps Capsulorhexis Utrata
- **41-4009** | Nagahara Chopper
- **11-6010** | Forceps Mcphearson Kelman Tying
- **41-5001** | Manipulator Lens Hook Sinskey Angled
- **20-7001** | Scissors Iris Straight Blunt
- **10-1106** | Forceps Bonn Straight Toothed Wrap 50x50 cm



Order Online: corzaeye.com

Order by Phone: 800-225-1195

For more information, visit corzaeye.com

Distributed by Corza Ophthalmology, (Katena Products Inc)

Corza Ophthalmology (Katena Products, Inc.), 6 Campus Drive, Suite 310, Parsippany, NJ 07054, USA
Phone: 973.989.1600, Toll free (USA only): 800.225.1195, Fax: 973.989.8175

All trademarks are owned by the company and its affiliates.

©2024 Corza Medical All Rights Reserved. COR-CATBLINK-062124-Rev0