

Dry Eye Solutions



Meibomian Gland Expression Maskin Probes Kit

Maskin® Meibomian Gland
Intraductal Starter Kit

07-6044 (3) 1mm and (2) 2mm Probes




Maskin® Meibomian Gland Intraductal Probes & Dilator

07-6040-1 (76µm OD) 1mm Probe, 10/Box

07-6040-2 (76µm OD) 2mm Probe, 10/Box

07-6040-4 (76µm OD) 4mm Probe, 10/Box



*Designed And Developed By MGD Innovations™
MASKIN® Is A Registered Trademark Of Steven L. Maskin, M.D.
**US Patent No. 9,510,844 & European Patent No. 2,222,355

Silicone Long-Term Plugs

SuperEagle®

The simple solution for dry eye

Supplied 2 individually packaged sterile plugs with inserter/dilator per twin pack

PART #	DIAMETER	LENGTH
31302	0.4 – 0.6mm (small)	1.4mm
31312	0.6 – 0.8mm (medium)	1.7mm
31322	0.8 and larger (large)	1.7mm



EaglePlug®

The original Tapered-Shaft™ plug

Supplied 2 individually packaged sterile plugs with inserter/dilator per twin pack

PART #	DIAMETER	LENGTH
30342	0.4mm	1.1mm
30352	0.5mm	1.2mm
30362	0.6mm	1.3mm
30372	0.7mm	1.4mm
30382	0.8mm	1.5mm



SuperFlex®

The comprehensive solution for dry eye

Supplied 2 individually packaged sterile plugs with inserter/dilator per twin pack

PART #	DIAMETER	LENGTH
31732	0.3mm	1.1mm
31742	0.4mm	1.1mm
31752	0.5mm	1.2mm
31762	0.6mm	1.3mm
31772	0.7mm	1.4mm
31782	0.8mm	1.5mm
31792	0.9mm	1.7mm
31702	1.0mm	1.8mm
31712	1.1mm	2.0mm
31722	1.2mm	2.0mm
31802	1.3mm	2.0mm



Punctal Occlusion Eagle Vision Punctum Plugs

Eagle FlexPlug®

The original flexible Tapered-Shaft™ plug

Supplied 2 individually packaged sterile plugs with inserter/dilator per twin pack

PART #	DIAMETER	LENGTH
30742	0.4mm	1.1mm
30752	0.5mm	1.2mm
30762	0.6mm	1.3mm
30772	0.7mm	1.4mm
30782	0.8mm	1.5mm
30792	0.9mm	1.6mm



Plug1™

The one-size solution for dry eye

Supplied 2 individually packaged sterile plugs with inserter/dilator per twin pack

PART #	DIAMETER	LENGTH
30012	0.5mm -0.8mm	1.8mm



Flow Controllers

Tapered Shaft™ Flow Controller

Designed to regulate lacrimal flow

Supplied 2 individually packaged sterile plugs with inserter/dilator per twin pack

PART #	DIAMETER	INSIDE DIAMETER	LENGTH
30652	0.5mm	0.2mm	1.4mm
30662	0.6mm	0.2mm	1.4mm
30672	0.7mm	0.2mm	1.7mm
30682	0.8mm	0.2mm	1.9mm



EaglePlug TearFlow™

Up to a 150% larger lumen than the Tapered Shaft™ Flow Controller

Supplied 2 individually packaged sterile plugs with inserter/dilator per twin pack

PART #	OUTSIDE DIAMETER	INSIDE DIAMETER	LENGTH
3005	0.5mm	0.2mm	1.2mm
3006	0.6mm	0.3mm	1.3mm
3007	0.7mm	0.3mm	1.4mm
3008	0.8mm	0.3mm	1.5mm
3009	0.9mm	0.4mm	1.6mm
3010	1.0mm	0.5mm	1.8mm



Temporary Plugs

DuraPlug®

Extended temporary inserts effective 60-180 days

Supplied 2 sterile inserts per pack; 10 packages per box

Total of 20 inserts

PART #	DIAMETER	LENGTH
0076	0.2mm	2.0mm
0077	0.3mm	2.0mm
0078	0.4mm	2.0mm
0079	0.5mm	2.0mm



Collagen

Effective 3-5 days; dissolves completely in 7-10 days

Supplied 6 inserts per sterile foam pack; 10 packs per box

Total of 60 inserts

PART #	DIAMETER	LENGTH
0025	Insert 0.2mm	2.0mm
0026	Insert 0.3mm	2.0mm
0027	Insert 0.4mm	2.0mm
0052	"Shorts" 0.2mm	1.6mm
0053	"Shorts" 0.3mm	1.6mm
0054	"Shorts" 0.4mm	1.6mm




Amniotic Membrane

AmbioDisk®

Clarify™ processed, dehydrated amniotic membrane allografts. The allografts serve as a safe covering or barrier applied to the ocular surface following repair and reconstruction procedure in the treatment of ocular diseases and/or abnormalities. The Clarify™ process retains the amniotic intermediate layer (IL) which is a vital part of the membrane barrier present between the amnion and chorion layers.

PART #	SIZE	MEMBRANE LAYERS
AD-5090	9mm	amnion & IL
AD-5012	12mm	amnion & IL
AD-5150	15mm	amnion & IL
AF-1150	15mm	amnion, IL & chorion

Clarify™ is a trademark of StimLabs LLC.



Tear Diagnostics

Color Bar™


Schirmer Tear Test

Supplied 2 standardized sterile strips per pouch

50 pouches per box; total of 100 strips

PART #

0039



Meibomian Gland Expressors

K5-5900 Arita Meibomian Gland Compressor

Broad, smoothly polished tips; designed to gently express meibum from inflamed or occluded meibomian glands.




08-01716 Maskin® Meibum Expressor

The special diamond knurling pattern on the stainless steel rollers prevents and enables the instrument to engage and roll against the smooth conjunctiva without abrading the mucosal surface.



AE-4521 Yoshitomi Meibomian Gland Compressor

This forceps has curved, polished jaws to aide in the expression of the meibomian glands in patients afflicted with MGD.



08-01720 Tarse Meibum Titanium Expressor Forceps

The jaws of this forceps form a "V", closer together at the bottom and wider at the top. When squeezed onto the eyelid, the forces applied deep within the gland are stronger and direct the oil upward and toward the gland orifice—this prevents oils from being pushed down further into the gland.





A series of horizontal dotted lines for writing notes, filling most of the page.

Order Online: corzaeye.com

Order by Phone: 800-225-1195

For more information, visit corzaeye.com



Corza Ophthalmology, (Katena Products, Inc.)
6 Campus Drive, Suite 310, Parsippany, NJ 07054, USA
Phone: 973.989.1600, Toll free (USA only): 800.225.1195, Fax: 973.989.8175
All trademarks are owned by the company and its affiliates.

©2022 Corza Medical

All Rights Reserved. KX-CAT-011819-Rev4